


REVISION HISTORY				
REV	DESCRIPTION	DATE	DR BY	CH BY
-	PROTOTYPE	12/2/2010	AFP	JO
A	UPDATED REV'S ON PARTS LIST	2/20/2015	AFP	JO
B	ADDED ENGRAVING TO COVER SLEEVE, ADDED SHEET 3, REVISED HARDWARE PART NUMBERS	5/7/2015	criley	JO
C	ADDED LENS SPECIFICATION TO SHEET 3	12/14/2015	criley	JO



OUTLINE



LENS SPECIFICATION:

MODEL NUMBER: 54-10-7
PART NUMBER: 1-21033
LENS TYPE: WATER IMMERSION
WAVELENGTH RANGE: 488-910nm
NUMERICAL APERTURE: 0.7
F#: 0.70
APERTURE: 10.0mm
EFL: 7.0mm
WORKING DISTANCE: 3.74mm WATER
FOV: 0.1
HOUSING: STAINLESS STEEL



Preliminary Drawing:
Drawing to be used for
Prototype and Quote
purposes ONLY

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS AND INCLUDE CHEMICALLY APPLIED OR PLATED FINISHES		PROPRIETARY AND CONFIDENTIAL THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF SPECIAL OPTICS. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF SPECIAL OPTICS IS PROHIBITED.		
SCALE:	4:1			
WEIGHT:	168.68 g			
DR BY:	criley			
CH BY:	JO			
PROJECT #:	E2741	ASSEMBLY (54-10-7 @ 488-910nm)		
DATE DR:	1/9/2014			
DO NOT SCALE DRAWING		SIZE B	1-21033 SHEET 3 OF 3	REV C